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general information: Whitepeaks Coprplex PMMA Discs can be machined with any milling system. If required, data must be entered into the CAM module. Milling tools required are burs made for PMMA milling.

storage: Protect from heat and direct sunlight.

indication: Whitepeaks Coprplex PMMA Discs are not certified for intraoral use and therefor must not be used in the patients mouth. The material is made for the dental casting technique.

finishing: To cut out the milled restarations from the Coprplex disc, use burs or discs for cutting acrylic materials.

sprues and runners: Please check the rules for your casting alloys or pressable ceramics to design sprues and runners for casting. Acrylic molds for casting should be covered with a thin layer of wax before embedding.

embedding: Use a ring free embedding system or coat the ring with a ceramic ring liner before filling with investment. Mix the investment material using vacuum and according to instructions for use of the manufacturer. Heat the mold conventionally or with a speed heating system up to 580°C, holding time 30 min. Afterwards heat the mould to the required end temperature according to the used alloy or press ceramic. Follow the manufacturers instructions for use of the investment material, the alloy and / or press ceramic used.

disposal: Whitepeaks Coprplex PMMA Discs are made from pure PMMA. Waste can be recycled with other acrylics or can be disposed of through domestic residual waste system. PMMA is insoluble in water and is no hazard to groundwater.