

Zirconia	Liquids					
	Copran Color <small>(dipping liquid)</small>	Copran Monolith Color Base <small>(dipping liquid)</small>	Copran Monolith Color PaintOn	CopraSupreme Color Base <small>(dipping liquid)</small>	CopraSupreme Color PaintOn	CopraSmile Color PaintOn
Copran Zri	✓	✗	✗	✗	✗	✗
Copran Zri Monolith	✗	✓	✓	✗	✗	✗
CopraSupreme	✗	✗	✗	✓	✓	✗
CopraSmile	✗	✗	✗	✗	✗	✓