




# SHERAPRINT-PLATE

ARTICULATION AID SET FOR PRINTED MODELS






# SHERAPRINT-PLATE

*patent pending*

Articulation aid set for printed models

- Articulate and remove models easily and exactly
- Precise fit of the model in its original position every time
- Stable hold and easy handling
- Suitable for any articulator
- Set consists of data set for direct impression of the base plate to the model, base plate, magnet, magnetic cup, adhesive disc and screw



SHERA. *The model.*

SHERA Werkstoff-Technologie GmbH

Espohlstr. 53 | 49448 Lemförde | Germany | Tel.: +49 (0) 5443-9933-0 | [www.shera.de](http://www.shera.de)